Glutamine Assay Kit-WST® (货号: G268)

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1. ER Stress Cooperates with Hypernutrition to Trigger TNF-Dependent Spontaneous HCC Development, *Cancer Cell*, **2014**, 26, 331-343
2. Rhythmic oscillations of the microRNA miR-96-5p play a neuroprotective role by indirectly regulating glutathione levels, *Nature Communications*, **2014**, 5, 3823
3. Phospholipid methylation controls Atg32-mediated mitophagy and Atg8 recycling, *The EMBO Journal*, **2015**, 34(21), 2703-2719
4. Administering xCT inhibitors based on circadian clock improves antitumor effects, *Cancer Research*, **2017**, 77(23), 6603-6613
5. Cold stress-induced ferroptosis involves the ASK1-p38 pathway, *EMBO reports*, 2017, 18(11), 2067-2078
6. Chloroplast DNA Replication Is Regulated by the Redox State Independently of Chloroplast Division in *Chlamydomonas reinhardtii*, *Plant Physiology*, **2013**, 161, 2102-2112
7. Oxidative stress inhibits healthy adipose expansion through suppression of SREBF1-mediated lipogenic pathway, *Diabetes*, **2018**, 67(6), 1113-1127
8. Endogenous reactive oxygen species cause astrocyte defects and neuronal dysfunctions in the hippocampus: a new model for aging brain, *Aging Cell*, **2017**, 16, 39-51
9. Ubiquinol-10 Supplementation Activates Mitochondria Functions to Decelerate Senescence in Senescence-Accelerated Mice, *Antioxidants & Redox Signaling*, **2014**, 20(16), 2606-2620
10. Skeletal muscle-specific eukaryotic translation initiation factor 2 α phosphorylation controls amino acid metabolism and fibroblast growth factor 21-mediated non-cell-autonomous energy metabolism, *The FASEB Journal*, **2016**, 30(2), 798-812
11. Riluzole exerts distinct antitumor effects from a metabotropic glutamate receptor 1-specific inhibitor on breast cancer cells, *Oncotarget*, **2017**, 8(27), 44639-44653

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WST® Solution	1 ml × 1瓶
Enzyme Solution	20 μl × 1瓶
Buffer Solution	11 ml × 2瓶
Dilution Buffer	10 ml × 1瓶

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1. Manganese superoxide dismutase polymorphism affects the oxidized low-density lipoprotein-induced apoptosis of macrophages and coronary artery disease, *European Heart Journal*, **2008**, *29*, 1267-1274
2. Extracellular Superoxide Dismutase Accelerates Endothelial Recovery and Inhibits In-Stent Restenosis in Stented Atherosclerotic Watanabe Heritable Hyperlipidemic Rabbit Aorta, *Journal of the American College of Cardiology*, **2007**, *50*, 2249-2253
3. An Abnormal Mitochondrial-Hypoxia Inducible Factor-1-Kv Channel Pathway Disrupts Oxygen Sensing and Triggers Pulmonary Arterial Hypertension in Fawn Hooded Rats: Similarities to Human Pulmonary Arterial Hypertension, *Circulation*, **2006**, *113*, 2630-2641
4. Deficient plastidic fatty acid synthesis triggers cell death by modulating mitochondrial reactive oxygen species, *Cell Research*, **2015**, *25*, 621-633
5. A Tau-Targeted Multifunctional Nanocomposite for Combinational Therapy of Alzheimer's Disease, *ACS Nano*, **2018**, *12*(2), 1321-1338
6. Endothelial Cell Palmitoylproteomic Identifies Novel Lipid-Modified Targets and Potential Substrates for Protein Acyl Transferases, *Circulation Research*, **2012**, *110*(10), 1336-44
7. Lung EC-SOD overexpression attenuates hypoxic induction of Egr-1 and chronic hypoxic pulmonary vascular remodeling, *American Journal of Respiratory and Critical Care Medicine*, **2008**, *295*, L422-L430
8. Inhaled Ethyl Nitrite Prevents Hyperoxia-impaired Postnatal Alveolar Development in Newborn Rats, *American Journal of Respiratory and Critical Care Medicine*, **2007**, *176*, 291-299
9. In vivo guiding nitrogen-doped carbon nanozyme for tumor catalytic therapy, *Nature Communications*, **2018**, *9*(1), 1440
10. A lytic polysaccharide monooxygenase-like protein functions in fungal copper import and meningitis, *Nature Chemical Biology*, **2020**, doi.org/10.1038/s41589-019-0437-9
11. Nestin regulates cellular redox homeostasis in lung cancer through the Keap1-Nrf2 feedback loop, *Nature Communications*, **2019**, *10*, 5043
12. Allelopathic interactions of linoleic acid and nitric oxide increase the competitive ability of *Microcystis aeruginosa*, *The ISME Journal*, **2017**, *11*, 1865-1876
13. Phase I Study of Copper-Binding Agent ATN-224 in Patients with Advanced Solid Tumors, *Clinical Cancer Research*, **2008**, *14*, 7526-7534
14. Chemoprevention of Carcinogenic Progression to Esophageal Adenocarcinoma by the Manganese Superoxide Dismutase Supplementation, *Clin. Cancer Res.*, **2007**, *13*, 5176-5182
15. Development of insulin resistance and obesity in mice overexpressing cellular glutathione peroxidase, *Proc.Natl.Acad.Sci.*, **2004**, *101*, 8852-8857

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1. Contradictory Effects of Mitochondria and Non-Mitochondria-Targeted Antioxidants on Hepatocarcinogenesis by Altering DNA Repair in Mice, *Hepatology*, 2018, 67(2), 623-635
2. Green monomeric photosensitizing fluorescent protein for photo-inducible protein inactivation and cell ablation, *BMC Biology*, 2018, 16, 50
3. Sonodynamic therapy inhibits palmitate-induced beta cell dysfunction via PINK1/Parkin-dependent mitophagy, *Cell Death & Disease*, 2019, 10, 457
4. Blue-Emitting Electron-Donor/Acceptor Dyads for Naked-Eye Fluorescence Detection of Singlet Oxygen, *ChemPhotoChem*, 2017, 1(7), 299-303
5. Physical Plasma Elicits Immunogenic Cancer Cell Death and Mitochondrial Singlet Oxygen, *IEEE Transactions on Plasma Science*, 2018, 2(2), 138-146

DNA Damage Quantification Kit-AP Site Counting- (货号: DK02)

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Standard ARP-DNA (60 µl/well) or sample DNA (60 µl/well) + DNA Binding Solution (100 µl/well)

[20 samples]

	1	2	3	4	5	6	7	8	9	10	11	12	
A		0 ARP DNA std.			sample 1		sample 5			sample 13			
B		2.5 ARP DNA std.						sample 6			sample 14		
C		5 ARP DNA std.			sample 2		sample 7			sample 15			
D		10 ARP DNA std.						sample 8			sample 16		
E		20 ARP DNA std.						sample 9			sample 17		
F		40 ARP DNA std.					sample 10			sample 18			
G		blank		sample 3			sample 11			sample 19			
H				sample 4			sample 12			sample 20			


blank: 60 µl TE Buffer + 100 µl DNA Binding Solution/well
Do not use column 1 and 12

↓
Washing (250 µl/well)
HRP-Streptavidin (150 µl/well)
Washing (250 µl/well)
Substrate (100 µl/well)

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Washing HRP-Streptavidin and Substrate

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1. Sirt3 Mediates Reduction of Oxidative Damage and Prevention of Age-Related Hearing Loss under Caloric Restriction, *Cell*, **2010**, *143*, 802-812
2. S-adenosylmethionine regulates apurinic/aprimidinic endonuclease 1 stability: implication in hepatocarcinogenesis, *Gastroenterology*, **2009**, *136*(3), 1025-36
3. Impaired Integrity of DNA After Recovery From Inflammation Causes Persistent Dysfunction of Colonic Smooth Muscle, *Gastroenterology*, **2011**, *141*, 1293-1301
4. DNA damage is a novel response to sublytic complement C5b-9-induced injury in podocytes, *The Journal of Clinical Investigation*, **2003**, *111*, 877-885
5. Symbiophagy as a cellular mechanism for coral bleaching, *Autophagy*, **2009**, *5*(2), 211-6
6. The prion protein is critical for DNA repair and cell survival after genotoxic stress, *Nucleic Acids Research*, **2015**, *43*(2), 904-916
7. APE1 deficiency promotes cellular senescence and premature aging features, *Nucleic Acids Research*, **2018**, *46*(11), 5664-5677
8. Telomeric epigenetic response mediated by Gadd45a regulates stem cell aging and lifespan, *EMBO Reports*, **2018**, *19*, e45494
9. Inhibiting translesion DNA synthesis as an approach to combat drug resistance to DNA damaging agents, *Oncotarget*, **2017**, *8*(25), 40804-40816

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
1. Cardiac Myocyte-Specific Expression of Inducible Nitric Oxide Synthase Protects Against Ischemia/Reperfusion Injury by Preventing Mitochondrial Permeability Transition, *Circulation*, **2008**, 118(19),1970-8
2. Genetic Deficiency of Glutathione S-Transferase P Increases Myocardial Sensitivity to Ischemia-Reperfusion Injury, *Circulation Research*, **2015**, 117(5), 437-49
3. Prussian Blue Nanoparticles as Multienzyme Mimetics and Reactive Oxygen Species Scavengers, *Journal of the American Chemical Society*, **2016**, 138(18), 5860-5
4. Highly bioactive zeolitic imidazolate framework-8-capped nanotherapeutics for efficient reversal of reperfusion-induced injury in ischemic stroke, *Science Advances*, **2020**, 6(12), eaay9751
5. Potentiating antibiotics in drug-resistant clinical isolates via stimuli-activated superoxide generation, *Science Advances*, **2017**, 3(10), e1701776
6. Detection and imaging of the free radical DNA in cells-Site-specific radical formation induced by Fenton chemistry and its repair in cellular DNA as seen by electron spin resonance, immuno-spin trapping and confocal microscopy, *Nucleic Acids Res.*, **2012**, 40(12), 5477-5486
7. Ultrasmall Cu₂-xS nanodots as photothermal-enhanced Fenton nanocatalysts for synergistic tumor therapy at NIR-II biowindow, *Biomaterials*, **2019**, 206, 101-114
8. Intelligent Nanocomposites with Intrinsic Blood-Brain-Barrier Crossing Ability Designed for Highly Specific MR Imaging and Sonodynamic Therapy of Glioblastoma, *Small*, **2020**, e1906985
9. Amino acid oxidation of the D1 and D2 proteins by oxygen radicals during photoinhibition of Photosystem II, *Proc. Natl. Acad. Sci.*, **2017**, 114(11), 2988-2993
10. Hydrogen sulfide cytoprotective signaling is endothelial nitric oxide synthase-nitric oxide dependent, *Proc. Natl. Acad. Sci.*, **2014**, 111(8), 3182-3187
11. Nonenzymatic displacement of chlorine and formation of free radicals upon the reaction of glutathione with PCB quinones, *Proc. Natl. Acad. Sci.*, **2009**, 106(24), 9725-9730
12. Overexpression of Extracellular Superoxide Dismutase Attenuates Heparanase Expression and Inhibits Breast Carcinoma Cell Growth and Invasion, *Cancer Research*, **2009**, 69(15), 6355-63
13. Manganese Superoxide Dismutase Modulates Hypoxia-Inducible Factor-1A Induction via Superoxide, *Cancer Research*, **2008**, 68(8), 2781-8
14. Iron and sulfur co-doped graphite carbon nitride (FeOy/S-g-C₃N₄) for activating peroxymonosulfate to enhance sulfamethoxazole degradation, *Chemical Engineering Journal*, **2020**, 382, 122836
15. Spontaneous DNA damage to the nuclear genome promotes senescence, redox imbalance and aging, *Redox Biology*, **2018**, 17, 259-273
16. Switch of Mitochondrial Superoxide Dismutase into a Prooxidant Peroxidase in Manganese-Deficient Cells and Mice, *Cell Chemical Biology*, **2018**, 25, 413-425
17. Photochemical reaction of tricresyl phosphate (TCP) in aqueous solution: Influencing factors and photolysis products, *Chemosphere*, **2020**, 241, 124971
18. Crossover between anti- and pro-oxidant activities of different manganese oxide nanoparticles and their biological implications, *Journal of Materials Chemistry B*, **2020**, DOI:10.1039/c9tb02524c

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名称: 5-tert-Butoxycarbonyl-5-methyl-1-pyrroline-N-oxide

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
- Highly Catalytic Niobium Carbide (MXene) Promotes Hematopoietic Recovery after Radiation by Free Radical Scavenging, *ACS Nano*, **2019**, 13(6), 6438-6454
- Peroxymonosulfate enhanced visible light photocatalytic degradation bisphenol A by single-atom dispersed Ag mesoporous g-C₃N₄ hybrid, *Applied Catalysis B: Environmental*, **2017**, 211, 79-88
- Synergetic Activation of Peroxymonosulfate by Co₃O₄ Modified g-C₃N₄ for Enhanced Degradation of Diclofenac Sodium under Visible Light Irradiation, *Applied Catalysis B: Environmental*, **2017**, 218, 810-818
- Tailored synthesis of active reduced graphene oxides from waste graphite: Structural defects and pollutant-dependent reactive radicals in aqueous organics decontamination, *Applied Catalysis B: Environmental*, **2018**, 229, 71-80
- Mitochondria-targeted TPP-MoS₂ with dual enzyme activity provides efficient neuroprotection through M1/M2 microglial polarization in an Alzheimer's disease model, *Biomaterials*, **2020**, 232, 119752
- Impact of humic acid on the degradation of levofloxacin by aqueous permanganate: Kinetics and mechanism, *Water Research*, **2017**, 123, 67-74
- Quinone group enhances the degradation of levofloxacin by aqueous permanganate: Kinetics and mechanism, *Water Research*, **2018**, 143, 109-116
- Enhanced Fenton-like degradation of pharmaceuticals over framework copper species in copper-doped mesoporous silica microspheres, *Chemical Engineering Journal*, **2015**, 274, 298-306
- MOF-templated synthesis of CoFe₂O₄ nanocrystals and its coupling with peroxymonosulfate for degradation of bisphenol A, *Chemical Engineering Journal*, **2018**, 353, 329-339
- Mechanism of Catalytic Ozonation in Fe₂O₃/Al₂O₃@SBA-15 Aqueous Suspension for Destruction of Ibuprofen, *Environ. Sci. Technol.*, **2015**, 49(3), 1690-1697
- Mechanism for enhanced degradation of clofibrac acid in aqueous bycatalytic ozonation over MnOx/SBA-15Qiangqiang, *Journal of Hazardous Materials*, **2015**, 286, 276-284