

连续分液器(Nichimate Stepper)

编号	型号	规格	设定容量(μL)	精度(%)	重复性(%)	数量	
2-5001-01	00-NSTP	分液器	-			1个	
2-5001-20	00-NS-STD	分液器支架					
2-5001-11	00-NS-V	0.05mL分液管(带吸头)	1/2/3/4/5	± 5	$< 8.0 \sim < 3.0$	1袋 (100个)	
2-5001-12	00-NS-S	0.5mL分液管	10/20/30/40/50	± 1.5	$< 1.4 \sim < 0.6$		
2-5001-13	00-NS-M	1.25mL分液管	25/50/75/100/125	± 1.2	$< 0.8 \sim < 0.5$		
2-5001-14	00-NS-N	2.5mL分液管	50/100/150/200/250	± 0.7	$< 0.8 \sim < 0.3$		
2-5001-15	00-NS-L	5.0mL分液管	100/200/300/400/500	± 0.5	$< 0.6 \sim < 0.3$		
2-5001-16	00-NS-X	12.5mL分液管	250/500/750/1000/1250				
2-5001-17	00-NS-Y	25mL分液管(带1个附件)	500/1000/1500/2000/2500	± 0.3	$< 0.6 \sim < 0.2$		1袋(50个)
2-5001-18	00-NS-Z	50mL分液管(带1个附件)	1000/2000/3000/4000/5000				1袋(25个)
61-0007-78	00-NS-SB	0.5mL分液管(已灭菌)	10/20/30/40/50	± 1.7	< 1.8	1袋 (100个) 独立 包装	
61-0007-76	00-NS-MB	1.25mL分液管(已灭菌)	25/50/75/100/125	± 1.6	< 2.1		
61-0007-77	00-NS-NB	2.5mL分液管(已灭菌)	50/100/150/200/250	± 1.2	< 1.7		
61-0007-75	00-NS-LB	5mL分液管(已灭菌)	100/200/300/400/500	± 2.0	< 2.4		
61-0007-79	00-NS-XB	12.5mL分液管(已灭菌)	250/500/750/1000/1250	± 1.1	< 1.7		